medskin solutions
DR. SUWELACK

MedSkin Solutions Dr. Suwelack AG
Josef-Suwelack-Strasse 2
48727 Billerbeck/Germany
phone: +49 2543 21820
Email: info@medskin-suwelack.com
www.medskin-suwelack.com
www.matriderm.com

medskin solutions
DR. SUWELACK

matriderm® plus+

BI-LAYER DERMAL MATRIX

تعليمات الاستخدام	**Használati utasítás**
Инструкции за употреба	**Istruzioni per l'uso**
Návod k použití	**Lietošanas pamācība**
Brugsanvisning	**Gebruiksaanwijzing**
Gebrauchsanweisung	**Bruksanvisning**
Οδηγίες χρήσης	**Instrukcja użycia**
Instructions For Use	**Instruções de utilização**
Instrucciones de uso	**Instrucțiuni de utilizare**
Kasutusjuhend	**Návod na použitie**
Käyttöohjeet	**Navodila za uporabo**
Notice d'utilisation	**Bruksanvisning**
Upute za uporabu	**Kullanım Talimatı**

AR 5–10	**تعليمات الاستخدام** اللغة العربية	HU 77–82	**Használati utasítás** magyar nyelvű
BG 11–16	**Инструкции за употреба** български език	IT 83–88	**Istruzioni per l'uso** lingua italiana
CS 17–22	**Návod k použití** čeština	LV 89–94	**Lietošanas pamācība** latviešu valoda
DA 23–28	**Brugsanvisning** dansk	NL 95–100	**Gebruiksaanwijzing** Nederlands
DE 29–34	**Gebrauchsanweisung** Deutsche Sprache	NO 101–106	**Bruksanvisning** norsk
EL 35–40	**Οδηγίες χρήσης** ελληνική γλώσσα	PL 107–112	**Instrukcja użycia** język polski
EN 41–46	**Instruction For Use** English Language	PT 113–118	**Instruções de utilização** português
ES 47–52	**Instrucciones de uso** español	RO 119–124	**Instrucțiuni de utilizare** limba română
ET 53–58	**Kasutusjuhend** eesti keel	SK 125–130	**Návod na použitie** slovenský jazyk
FI 59–64	**Käyttöohjeet** suomi	SL 131–136	**Navodila za uporabo** slovenski jezik
FR 65–70	**Notice d'utilisation** français	SV 137–142	**Bruksanvisning** svenska
HR 71–76	**Upute za uporabu** hrvatski jezik	TR 143–148	**Kullanım Talimatı** Türkçe

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matriderm® pluS+
BI-LAYER DERMAL MATRIX

Important Usage and Safety Information for MatriDerm pluS+

Reporting Serious Incidents: Any serious incident involving MatriDerm pluS+ must be reported to MedSkin Solutions Dr. Suwelack AG and the competent authority of the Member State where the user and / or patient is located.

Animal Origin: MatriDerm pluS+ contains animal-derived tissues. The risk of prion contamination is theoretically possible but highly unlikely and minimal.

Single-use only: MatriDerm pluS+ is sterilized using gamma irradiation and is intended for single use only. Before use, ensure the sterile packaging is undamaged and unopened; do not use if compromised. Once opened, the product cannot be stored or reused. Discard any unused pieces. Re-using the product poses risks of transmitting fluids or tissue particles containing pathogens from previous patients. Do not resterilize MatriDerm pluS+, because this can affect the physicochemical properties of the matrix.

Expiration Date: Do not use MatriDerm pluS+ after the expiration date stated on the packaging.

Patients at increased risk: Exercise caution with patients undergoing radiotherapy and/or cytotoxic treatments (e.g., chemotherapy), as they may experience delayed wound healing. Patients with untreated metabolic disorders (e.g., diabetes, osteomalacia, thyroid disease), long-term corticosteroid use, autoimmune conditions, or heavy smoking habits are also at risk of delayed wound healing.

Incompatible Substances: Do not let MatriDerm pluS+ come into direct contact with iodine-containing antiseptics, enzymatic debridement agents, collagenases, or caustic agents that alter proteins.

If these substances are used, follow their respective Instructions For Use (IFUs) and ensure all residuals are thoroughly rinsed off before applying MatriDerm pluS+.

Tissue Specificity: Do not apply on necrotic tissue.

No Pregnancy Data: No clinical data are available concerning use during pregnancy or breastfeeding.

Disposal: Unused product or parts thereof may be disposed of in accordance with institutional procedures for handling biological waste.

Storage Conditions and Shelf Life

Storage: MatriDerm pluS+ should be kept in a dry place at room temperature (do not store above 25°).

Instructions For Use

Shelf Life: The shelf life of MatriDerm pluS+ (unopened in original packaging) is 36 months and it should not be used after the stated expiry date. The expiry date is printed on the packaging.

Intended Use

MatriDerm pluS+ consists of two layers: The first layer consists of a silicone grid. The second layer is a three-dimensional fenestrated matrix composed of native structurally intact collagen fibrils and elastin. For the reconstruction of dermal tissue MatriDerm pluS+ serves as a dermal substitute and may be used in combination with a split-thickness skin graft (STSG) in a second stage, after prior removal of the silicone grid. The aim of treatment is to construct a neo dermis, to improve the quality of reconstructed tissue, reduce scarring, prevent wound contraction and restore functionality.

Intended Users: MatriDerm pluS+ is to be used by qualified healthcare professionals.

Intended User Environment: MatriDerm pluS+ is to be used in a sterile environment in an operating room.

Patient Population: MatriDerm pluS+ can be applied to pediatric and adult patients in need of dermis reconstruction except for pregnant and/or lactating and breast-feeding women.

Indications for Use

MatriDerm pluS+ is indicated for the treatment of wounds including:

- Partial thickness and full-thickness wounds
- Graft-requiring poorly healing (e.g. chronic) wounds (such as pressure ulcers, venous ulcers, diabetic ulcers, chronic ulcers)
- Surgical wounds (dermatological diseases, donor sites/grafts, post Mohs surgery, post-laser surgery, podiatric, wound dehiscence)
- 2nd degree and 3rd degree burns
- Trauma wounds (abrasions, lacerations, and skin tears)
- Draining wounds
- Mucosal defects

Contraindications

Hypersensitivity: Contraindicated in patients with known hypersensitivity to bovine collagen, elastin, or silicone.

Infected wounds: Do not use on clinically diagnosed infected wounds.

Awareness of Risks: Inform patients about potential risks, including infection, failure to incorporate into the wound bed, partial or complete graft loss,

seromas, hematomas, shearing or mechanical loss of the product, delayed healing and hyperpigmentation. The patient should be made aware of these risks and others associated with this kind of surgery.

Intolerance to Collagen: Despite the low antigenic potential of collagen described in the literature, an intolerance reaction cannot be entirely excluded.

Collagenases and enzymatic agents: Simultaneous use of collagenases, such as those used for enzymatic wound cleansing, can accelerate the degradation of MatriDerm pluS+.

Patient Implant Card

An implant card and a patient sticker with device specific information (e.g. LOT number) will be provided with the product packaging. The implant card must be completed by the healthcare Institution/provider by filling out the following patient specific details:

- Name of the patient or patient ID
- Date of implantation
- Name and address of the healthcare institution/provider.

Place one patient sticker into the according field in the implant card.

Clinical Benefit

MatriDerm pluS+ enables closure of partial and full thickness wounds in a two-step procedure and results in a neo dermis with a good quality of the grafted skin. Furthermore, clinical trials on the treatment of full thickness skin defects with MatriDerm Dermal Matrix have shown that the elasticity of the regenerated skin is similar to that of the perilesional skin. MatriDerm pluS+ can be used over exposed tendons and bones (non-critical size defects) to cover the defect with placement of a STSG.

Summary of Safety and Clinical performance

For further information on the safety and clinical performance of this device please see the summary of safety and clinical performance (SSCP) in the EUDAMED database:
https://ec.europa.eu/tools/eudamed/
Basic-UDI-DI: 42505690114V

Device Composition and Origin

Device Composition: MatriDerm pluS+ consists of approximately 96.7% bovine collagen and 2.3% bovine elastin plus a medical grade silicone grid.

MatriDerm pluS+ is a natural product; dimensions are approximate.

Before Use

Inspect Packaging: Verify that the sterile pouch is intact and not past its expiry date; discard if the package is damaged or previously opened.

Single application only: Do not attempt to resterilize.

Method of Application

1. Preparation

- **Ensure Aseptic Technique:** Wash hands thoroughly and wear sterile gloves. Use sterile instruments and materials.
- **Clean the Wound:** Gently clean the wound with a sterile physiological saline solution to remove any debris or contaminants.
- **Debride the Wound:** Perform radical debridement to remove all non-viable tissue. This step is crucial for promoting optimal healing.
- **Achieve Hemostasis:** Stop any active bleeding by applying gentle pressure or using appropriate hemostatic agents.
- **Prepare the Wound Bed:** Ensure the wound bed is non-bleeding but still moist. This promotes better integration of the collagen matrix.

2. Placement:
Ensure the collagen-elastin layer contacts the wound directly, with the silicone layer facing upwards.

3. Cut to Size: Trim MatriDerm pluS+ to match the wound dimensions, with a slightly larger overlap (1-2 mm) on the wound bed edges.

4. Moistening: Moisten the collagen matrix and the non-bleeding wound with a physiological saline solution at room temperature. Ensure hemostasis without damaging the wound bed.

5. Secure Collagen Matrix: Smooth out any air bubbles with gentle, even pressure to avoid collapse and ensure proper vascularization; firmly secure MatriDerm pluS+ to the wound using surgical staples or equivalent. MatriDerm pluS+ is compatible with all negative pressure wound therapy (NPWT) devices.

Note: Steps 3 & 4 can also be reversed, i.e. the matrix can be moistened first and then be cut to size.

Post-Application Care

Immediately: After Application
Apply Secondary Dressing: Choose a secondary dressing suitable for the wound's location, size, and depth. This dressing will help maintain MatriDerm pluS+ in place, ensure a moist wound environment and protect the area from external contaminants.

Day 1 to 3: Initial Monitoring
Inspect for Complications: Check the

wound for signs of infection, failure to incorporate into the wound bed, partial or complete graft loss, hematoma, seromas, shearing or mechanical loss of the product, delayed healing and hyperpigmentation, or excessive fluid buildup.

Drainage: If fluid buildup occurs, create small openings in the collagen matrix for drainage.

Day 3 to 5: First Dressing Change

Leave Silicone Layer in Place Initially: Leave silicone layer in place during the first dressing change to avoid disturbing the wound bed.

Frequency of Changes: Adjust how often you change the dressing based on the amount of fluid coming from the wound and the wound's condition.

Day 5 to 12: Silicone Layer Removal

Assess Vascularization: Between 5 to 12 - days post-application, observe the wound for signs of vascularization (new blood vessel formation).

Remove Silicone Layer Gently: When vascularization is evident, gently remove the silicone layer.

Don't Remove the Collagen-Elastin Layer: MatriDerm pluS+ is single-use and bio-degrades over 4-6 weeks. No replacement is needed.

Reduce Granulation tissue: Excess granulation tissue that has extended beyond the matrix should be reduced to produce a suitable wound bed for the graft. Hemostasis should be achieved.

Apply STSG: after harvesting a 0,20 mm from the donor site apply the Split Thickness Skin Graft (meshed or unmeshed) on the surface of the wound and secure it using surgical staples or equivalent.

Apply Secondary Dressing: To ensure a moist wound environment, apply a secondary dressing according to hospital protocol, with an interface of a non-adherent wound dressing, and monitor closely.

Up to 6 weeks: Monitor Healing

Inspect Regularly: Monitor for infection, healing, and complications. Adjust secondary dressings to manage exudate and protect the wound.

Explanation of Symbols on Packaging

Use-by date		Manufacturer	
Upper limit of temperature		Contains biological material of animal origin	
Keep dry		Caution	
Keep away from sunlight		Medical device	
Sterilized using irradiation		Unique device identifier	
Do not re-use		Double sterile barrier system	
Do not resterilize		Patient identification	
Do not use if package is damaged and consult instructions for use		Patient information website	
CE-sign with identification number of the Notified Body		Date	
Batch Code		Health care centre or doctor	
Catalogue number		Quantity	
Consult electronic instructions for use		Thickness of the product	
Date of manufacture		Humidity limitation (upper limit)	

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